



Introducing VerraGlo™ Oil Oxidation Test for Fats & Oils

Simplified, Rapid Peroxide Value Testing

Do oxidized oils give your food products off odors and flavors? Are you specifying peroxide values in the oils you are purchasing for processing? Does difficult iodometric titration measurement make verification testing difficult?


Charm Sciences, Inc. introduces a revolutionary **Peroxide Value by Light Emitting Diode (PV-LED)** assay, read on the handheld novaLUM II-X luminometer. *This assay will take your oil quality testing into the future.*

Benefits Include:

- 2-minute, simple procedure - no formal training required
- Ready-to-use reagent can be performed in non-laboratory setting
- No toxic chemicals or hazardous waste
- No glassware or clean up - disposable reagents
- Improved sensitivity, precision, and readability in the 0.3-20 PV range



How it works



Scan the QR code to speak to our Charm experts!

www.charm.com

T +1.978.687.9200 | F +1.978.687.9216 | E info@charm.com | 659 Andover Street, Lawrence, MA 01843

© 2023 Charm Sciences, Inc. | Charm is a registered trademark of Charm Sciences, Inc. | MRK-1496-001 Mar-2023